



# CD88 Monoclonal Antibody

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>         | YP-mAb-14093                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Isotype</b>            | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Reactivity</b>         | Human;Rat;Mouse;                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Applications</b>       | WB                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Gene Name</b>          | C5AR1                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Protein Name</b>       | C5a anaphylatoxin chemotactic receptor                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Immunogen</b>          | The antiserum was produced against synthesized peptide derived from the N-terminal region of human C5AR1. AA range:1-50                                                                                                                                                                                                                                                                                                                                                   |
| <b>Specificity</b>        | CD88 Monoclonal Antibody detects endogenous levels of CD88 protein.                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Formulation</b>        | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Source</b>             | Monoclonal, Mouse,IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purification</b>       | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                                                                                                                                      |
| <b>Dilution</b>           | WB 1:500-1:2000                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Concentration</b>      | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Purity</b>             | ≥90%                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Storage Stability</b>  | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Synonyms</b>           | C5AR1; C5AR; C5R1; C5a anaphylatoxin chemotactic receptor; C5a-R; C5aR; CD88                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Observed Band</b>      | 39kD                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Cell Pathway</b>       | Cell membrane ; Multi-pass membrane protein . Cytoplasmic vesicle . Phosphorylated C5aR colocalizes with ARRB1 and ARRB2 in cytoplasmic vesicles. .                                                                                                                                                                                                                                                                                                                       |
| <b>Tissue Specificity</b> | Brain,                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Function</b>           | function:Receptor for the chemotactic and inflammatory peptide anaphylatoxin C5a. This receptor stimulates chemotaxis, granule enzyme release and superoxide anion production.,PTM:Sulfation plays a critical role in the association of the receptor with C5a, but no significant role in the ability of the receptor to transduce a signal and mobilize calcium in response to a small peptide agonist.,similarity:Belongs to the G-protein coupled receptor 1 family., |
| <b>Background</b>         | function:Receptor for the chemotactic and inflammatory peptide anaphylatoxin C5a. This receptor stimulates chemotaxis, granule enzyme release and superoxide anion production.,PTM:Sulfation plays a critical role in the association of the receptor with C5a, but no significant role in the ability of the receptor to transduce a signal and mobilize calcium in response to a small peptide agonist.,similarity:Belongs to the G-protein coupled receptor 1 family., |




**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**



Western Blot analysis of various cells using CD88 Monoclonal Antibody