




ADCK5 Monoclonal Antibody

Catalog No	YP-mAb-14646
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ADCK5
Protein Name	Uncharacterized aarF domain-containing protein kinase 5
Immunogen	The antiserum was produced against synthesized peptide derived from human ADCK5. AA range:151-200
Specificity	ADCK5 Monoclonal Antibody detects endogenous levels of ADCK5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADCK5; Uncharacterized aarF domain-containing protein kinase 5
Observed Band	66kD
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	Brain,Colon,Ovary,Prostate,Small intestine,
Function	function:The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).,similarity:Belongs to the protein kinase superfamily. ADCK protein kinase family.,similarity:Contains 1 protein kinase domain.,
Background	function:The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).,similarity:Belongs to the protein kinase superfamily. ADCK protein kinase family.,similarity:Contains 1 protein kinase domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ADCK5 Monoclonal Antibody