




NARFL Monoclonal Antibody

Catalog No	YP-mAb-14873
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NARFL
Protein Name	Cytosolic Fe-S cluster assembly factor NARFL
Immunogen	The antiserum was produced against synthesized peptide derived from human NARFL. AA range:387-436
Specificity	NARFL Monoclonal Antibody detects endogenous levels of NARFL protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NARFL; PRN; Cytosolic Fe-S cluster assembly factor NARFL; Iron-only hydrogenase-like protein 1; IOP1; Nuclear prelamin A recognition factor-like protein; Protein related to Narf
Observed Band	50kD
Cell Pathway	CIA complex,
Tissue Specificity	Widely expressed.
Function	function:Regulates HIF1A expression.,similarity:Belongs to the NARF family.,tissue specificity:Widely expressed.,
Background	function:Regulates HIF1A expression.,similarity:Belongs to the NARF family.,tissue specificity:Widely expressed.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using NARFL Monoclonal Antibody