



Dlx-5 Monoclonal Antibody

Catalog No	YP-mAb-15754
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DLX5
Protein Name	Homeobox protein DLX-5
Immunogen	The antiserum was produced against synthesized peptide derived from human DLX5. AA range:91-140
Specificity	Dlx-5 Monoclonal Antibody detects endogenous levels of Dlx-5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DLX5; Homeobox protein DLX-5
Observed Band	32kD
Cell Pathway	Nucleus .
Tissue Specificity	Embryo,Lung,Osteoblast,Placenta,
Function	developmental stage:First expressed in embryos at 8.5-9 days in facial and branchial arch mesenchyme, otic vesicles and frontonasal ectoderm around olfactory placodes, a day later expression is seen in the developing forebrain in primordia of the ganglionic eminence and ventral diencephalic regions. In day 12.5 embryos, expressed in the brain and bones, and also in all skeletal structures of midgestation embryos after the first cartilage formation. Expression remains unaltered in both brain and skeleton in day 15 embryos and slowly decreases in day 17 embryos.,similarity:Belongs to the distal-less homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,subunit:Interacts with G22P1 (Ku70).,
Background	This gene encodes a member of a homeobox transcription factor gene family similar to the Drosophila distal-less gene. The encoded protein may play a role in bone development and fracture healing. Mutation in this gene, which is located in a tail-to-tail configuration with another member of the family on the long arm of




chromosome 7, may be associated with split-hand/split-foot malformation. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using Dlx-5 Monoclonal Antibody