



HoxA5 Monoclonal Antibody

Catalog No	YP-mAb-15761
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HOXA5
Protein Name	Homeobox protein Hox-A5
Immunogen	Synthesized peptide derived from the N-terminal region of human HoxA5.
Specificity	HoxA5 Monoclonal Antibody detects endogenous levels of HoxA5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HOXA5; HOX1C; Homeobox protein Hox-A5; Homeobox protein Hox-1C
Observed Band	40kD
Cell Pathway	Nucleus.
Tissue Specificity	Ovary,
Function	developmental stage:Expressed during embryogenesis and in adult kidney.,function:Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. Also binds to its own promoter. Binds specifically to the motif 5'-CYNATTA[TG]Y-3'.,similarity:Belongs to the Antp homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,
Background	In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. Methylation of this gene may result in the loss of its expression and, since the encoded protein upregulates the tumor suppressor p53, this protein may play an important role in tumorigenesis. [provided by RefSeq, Jul 2008],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using HoxA5 Monoclonal Antibody