



PRX1 Monoclonal Antibody

Catalog No	YP-mAb-15800
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRRX1
Protein Name	Paired mesoderm homeobox protein 1
Immunogen	Synthesized peptide derived from PRX1 . at AA range: 30-110
Specificity	PRX1 Monoclonal Antibody detects endogenous levels of PRX1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRRX1; PMX1; Paired mesoderm homeobox protein 1; Homeobox protein PHOX1; Paired-related homeobox protein 1; PRX-1
Observed Band	27kD
Cell Pathway	Nucleus.
Tissue Specificity	
Function	function:Acts as a transcriptional regulator of muscle creatine kinase (MCK) and so has a role in the establishment of diverse mesodermal muscle types. The protein binds to an A/T-rich element in the muscle creatine enhancer.,similarity:Belongs to the paired homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,similarity:Contains 1 OAR domain.,
Background	The DNA-associated protein encoded by this gene is a member of the paired family of homeobox proteins localized to the nucleus. The protein functions as a transcription co-activator, enhancing the DNA-binding activity of serum response factor, a protein required for the induction of genes by growth and differentiation factors. The protein regulates muscle creatine kinase, indicating a role in the establishment of diverse mesodermal muscle types. Alternative splicing yields two isoforms that differ in abundance and expression patterns. [provided by RefSeq, Jul 2008],



matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Fig. 4: The regulatory mechanism of miR-1307-mediated inhibition of **PRRX1**.



Western Blot analysis of various cells using PRRX1 Monoclonal Antibody