



GBP1 Monoclonal Antibody


Catalog No	YP-mAb-16159
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GBP1
Protein Name	Interferon-induced guanylate-binding protein 1
Immunogen	The antiserum was produced against synthesized peptide derived from human GBP1. AA range:71-120
Specificity	GBP1 Monoclonal Antibody detects endogenous levels of GBP1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GBP1; Interferon-induced guanylate-binding protein 1; GTP-binding protein 1; GBP-1; HuGBP-1; Guanine nucleotide-binding protein 1
Observed Band	68kD
Cell Pathway	Cytoplasm . Golgi apparatus membrane; Lipid-anchor ; Cytoplasmic side. Cell membrane . Secreted . Cytoplasmic vesicle . Secreted from endothelial cells in the cerebrospinal fluid, upon bacterial challenge and independently of IFNG induction. Golgi membrane localization requires isoprenylation and the presence of another IFNG-induced factor.
Tissue Specificity	Placenta,Spleen,Uterus,
Function	function:Binds GTP, GDP and GMP.,induction:By interferon gamma during macrophage activation.,similarity:Belongs to the GBP family.,
Background	Guanylate binding protein expression is induced by interferon. Guanylate binding proteins are characterized by their ability to specifically bind guanine nucleotides (GMP, GDP, and GTP) and are distinguished from the GTP-binding proteins by the presence of 2 binding motifs rather than 3. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using GBP1 Monoclonal Antibody