




GRTP1 Monoclonal Antibody

Catalog No	YP-mAb-16166
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GRTP1
Protein Name	Growth hormone-regulated TBC protein 1
Immunogen	The antiserum was produced against synthesized peptide derived from human GRTP1. AA range:71-120
Specificity	GRTP1 Monoclonal Antibody detects endogenous levels of GRTP1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GRTP1; TBC1D6; Growth hormone-regulated TBC protein 1; TBC1 domain family member 6
Observed Band	35kD
Cell Pathway	intracellular,endomembrane system,
Tissue Specificity	Brain,Cervix carcinoma,Kidney epithelium,Ovary,Pancreas,Uterus,
Function	function:May act as a GTPase-activating protein for Rab family protein(s).,similarity:Contains 1 Rab-GAP TBC domain.,
Background	function:May act as a GTPase-activating protein for Rab family protein(s).,similarity:Contains 1 Rab-GAP TBC domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GRTP1 Monoclonal Antibody