



# Rab11-FIP4 Monoclonal Antibody


|                           |                                                                                                                                                                                                                                                                                                                                                                         |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>         | YP-mAb-16213                                                                                                                                                                                                                                                                                                                                                            |
| <b>Isotype</b>            | IgG                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Reactivity</b>         | Human;Mouse                                                                                                                                                                                                                                                                                                                                                             |
| <b>Applications</b>       | WB                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Gene Name</b>          | RAB11FIP4                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Name</b>       | Rab11 family-interacting protein 4                                                                                                                                                                                                                                                                                                                                      |
| <b>Immunogen</b>          | The antiserum was produced against synthesized peptide derived from human RAB11FIP4. AA range:452-501                                                                                                                                                                                                                                                                   |
| <b>Specificity</b>        | Rab11-FIP4 Monoclonal Antibody detects endogenous levels of Rab11-FIP4 protein.                                                                                                                                                                                                                                                                                         |
| <b>Formulation</b>        | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                 |
| <b>Source</b>             | Monoclonal, Mouse,IgG                                                                                                                                                                                                                                                                                                                                                   |
| <b>Purification</b>       | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                                    |
| <b>Dilution</b>           | WB 1:500-1:2000                                                                                                                                                                                                                                                                                                                                                         |
| <b>Concentration</b>      | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Purity</b>             | ≥90%                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Storage Stability</b>  | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                            |
| <b>Synonyms</b>           | RAB11FIP4; ARFO2; KIAA1821; Rab11 family-interacting protein 4; FIP4-Rab11; Rab11-FIP4; Arfophilin-2                                                                                                                                                                                                                                                                    |
| <b>Observed Band</b>      | 95kD                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Cell Pathway</b>       | Endosome . Cytoplasm, cytoskeleton, spindle . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Recycling endosome membrane ; Peripheral membrane protein. Cleavage furrow . Midbody . Cytoplasmic vesicle . Recruited to the cleavage furrow and the midbody during cytokinesis. .                                                                  |
| <b>Tissue Specificity</b> | Present at high level in testis (at protein level). Weakly expressed in other tissues.                                                                                                                                                                                                                                                                                  |
| <b>Function</b>           | function:A Rab11 effector protein acting in the endosomal trafficking.,similarity:Contains 1 EF-hand domain.,subcellular location:Colocalizes with RAB11A in the endocytic recycling compartment (ERC).,subunit:Homooligomer. Forms an heterooligomeric complex with RAB11FIP2, RAB11FIP3 and RAB11FIP5. Interacts with RAB11A that has been activated by GTP binding., |
| <b>Background</b>         | Proteins of the large Rab GTPase family (see RAB1A; MIM 179508) have regulatory roles in the formation, targeting, and fusion of intracellular transport vesicles. RAB11FIP4 is one of many proteins that interact with and regulate Rab GTPases (Hales et al., 2001 [PubMed 11495908]).[supplied by OMIM, Apr 2008],                                                   |

**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**

Western Blot analysis of various cells using Rab11-FIP4 Monoclonal Antibody