




RHOBTB3 Monoclonal Antibody

Catalog No	YP-mAb-16246
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RHOBTB3
Protein Name	Rho-related BTB domain-containing protein 3
Immunogen	The antiserum was produced against synthesized peptide derived from human RHOBTB3. AA range:335-384
Specificity	RHOBTB3 Monoclonal Antibody detects endogenous levels of RHOBTB3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RHOBTB3; KIAA0878; Rho-related BTB domain-containing protein 3
Observed Band	70kD
Cell Pathway	Golgi apparatus .
Tissue Specificity	Ubiquitous. Highly expressed in neural and cardiac tissues, pancreas, placenta and testis.
Function	similarity:Contains 2 BTB (POZ) domains.,tissue specificity:Ubiquitous. Highly expressed in neural and cardiac tissues, pancreas, placenta and testis.,
Background	RHOBTB3 is a member of the evolutionarily conserved RHOBTB subfamily of Rho GTPases. For background information on RHOBTBs, see RHOBTB1 (MIM 607351).[supplied by OMIM, Apr 2004],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RHOBTB3 Monoclonal Antibody