

SFTPC mouse mAb

货号	YP-mAb-17991
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	SFTPC
基因名称	SFTPC SFTP2
蛋白名称	Pulmonary surfactant-associated protein C (SP-C) (Pulmonary surfactant-associated proteolipid SPL(Val)) (SP5)
免疫原	Synthesized peptide derived from human SFTPC
特异性	This antibody detects endogenous levels of SFTPC at Human,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Monoclonal, mouse,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	22kD
功能	Pulmonary surfactant associated proteins promote alveolar stability by lowering the surface tension at the air-liquid interface in the peripheral air spaces.
细胞定位	Secreted, extracellular space, surface film.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images