



PRMT3 Rabbit pAb

Catalog No	YP-Ab-18878
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	PRMT3 HRMT1L3
Protein Name	Protein arginine N-methyltransferase 3 (Heterogeneous nuclear ribonucleoprotein methyltransferase-like protein 3)
Immunogen	Synthesized peptide derived from human PRMT3
Specificity	This antibody detects endogenous levels of PRMT3 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	58kD
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	Protein-arginine N-methyltransferase that catalyzes both the monomethylation and asymmetric dimethylation of the guanidino nitrogens of arginine residues in target proteins, and therefore falls into the group of type I methyltransferases (Probable). May regulate retinoic acid synthesis and signaling by inhibiting ALDH1A1 retinal dehydrogenase activity .
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images