



SAMD9L Rabbit pAb

Catalog No	YP-Ab-18882
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	SAMD9L C7orf6 DRIF2 KIAA2005 UEF
Protein Name	Sterile alpha motif domain-containing protein 9-like (SAM domain-containing protein 9-like)
Immunogen	Synthesized peptide derived from human SAMD9L
Specificity	This antibody detects endogenous levels of SAMD9L at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	174kD
Cell Pathway	Early endosome .
Tissue Specificity	Widely expressed in adult and fetal tissues. Variable expression in tumors. Down-regulated in breast cancer.
Function	May be involved in endosome fusion. Mediates down-regulation of growth factor signaling via internalization of growth factor receptors.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images