



SLC27A3 Rabbit pAb

Catalog No	YP-Ab-18883
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	
Protein Name	
Immunogen	Synthesized peptide derived from human SLC27A3
Specificity	This antibody detects endogenous levels of SLC27A3 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	75kD
Cell Pathway	Mitochondrion membrane ; Single-pass membrane protein .
Tissue Specificity	Expressed in bronchial and bronchiolar epithelial cells (at protein level).
Function	Has acyl-CoA ligase activity for long-chain and very-long-chain fatty acids . Does not exhibit fatty acid transport activity (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images