



P2RX7 rabbit pAb

Catalog No	YP-Ab-17274
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	IHC,WB
Gene Name	P2RX7
Protein Name	P2X purinoceptor 7 (P2X7) (ATP receptor) (P2Z receptor) (Purinergic receptor)
Immunogen	Synthesized peptide derived from human C-ternal P2RX7
Specificity	This antibody detects endogenous levels of P2RX7 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit,polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	P2X purinoceptor 7 (P2X7) (ATP receptor) (P2Z receptor) (Purinergic receptor)
Observed Band	65kD
Calculated Molecular Weight	33kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed with highest levels in brain and immune tissues.
Function	Receptor for ATP that acts as a ligand-gated ion channel. Responsible for ATP-dependent lysis of macrophages through the formation of membrane pores permeable to large molecules. Could function in both fast synaptic transmission and the ATP-mediated lysis of antigen-presenting cells. In the absence of its natural ligand, ATP, functions as a scavenger receptor in the recognition and engulfment of apoptotic cells .
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images