



FLAP rabbit pAb

货号	YP-Ab-17928
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	ALOX5AP
基因名称	ALOX5AP FLAP
蛋白名称	Arachidonate 5-lipoxygenase-activating protein (FLAP) (MK-886-binding protein)
免疫原	Synthesized peptide derived from human FLAP
特异性	This antibody detects endogenous levels of FLAP at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	18kD
功能	Required for leukotriene biosynthesis by ALOX5 (5-lipoxygenase). Anchors ALOX5 to the membrane. Binds arachidonic acid, and could play an essential role in the transfer of arachidonic acid to ALOX5. Binds to MK-886, a compound that blocks the biosynthesis of leukotrienes.
细胞定位	Nucleus membrane; Multi-pass membrane protein. Endoplasmic reticulum membrane; Multi-pass membrane protein.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

