




ZDH11 mouse mAb

Catalog No	YP-mAb-08628
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZDHHC11 ZNF399
Protein Name	ZDH11
Immunogen	Synthesized peptide derived from human ZDH11 AA range: 17-67
Specificity	This antibody detects endogenous levels of ZDH11 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1 : 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	45kD
Observed Band	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in testis.
Function	catalytic activity:Palmitoyl-CoA + protein-cysteine = S-palmitoyl protein + CoA.;domain:The DHHC domain is required for palmitoyltransferase activity.;similarity:Belongs to the DHHC palmitoyltransferase family.;similarity:Contains 1 DHHC-type zinc finger.;
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZDH11 mouse mAb