



ZDHC1 mouse mAb

Catalog No	YP-mAb-10011
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ZDHHC1 C16orf1 ZNF377
Protein Name	Probable palmitoyltransferase ZDHHC1 (EC 2.3.1.-) (DHHC domain-containing cysteine-rich protein 1) (Zinc finger DHHC domain-containing protein 1) (DHHC-1) (Zinc finger protein 377)
Immunogen	Synthesized peptide derived from human ZDHC1 AA range: 81-131
Specificity	This antibody detects endogenous levels of human ZDHC1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	53kD
Observed Band	
Cell Pathway	Endosome membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein . Golgi apparatus .
Tissue Specificity	Widely expressed with significant expression in heart, brain, placenta, lung, liver, kidney, testis, thymus and small intestine (PubMed:16647879). Expressed at lower levels in adult pancreas and lung (PubMed:10395086).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images