



ACYP1 Rabbit pAb

Catalog No	YP-Ab-18922
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	ACYP1 ACYPE
Protein Name	Acylphosphatase-1 (Acylphosphatase, erythrocyte isozyme) (Acylphosphatase, organ-common type isozyme) (Acylphosphate phosphohydrolase 1)
Immunogen	Synthesized peptide derived from human ACYP1
Specificity	This antibody detects endogenous levels of ACYP1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	11kD
Cell Pathway	
Tissue Specificity	Organ-common type isozyme is found in many different tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images