



RS16 rabbit pAb

Catalog No	YP-Ab-09584
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPS16
Protein Name	40S ribosomal protein S16
Immunogen	Synthesized peptide derived from human RS16 AA range: 39-89
Specificity	This antibody detects endogenous levels of human RS16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Calculated Molecular Weight	16kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

