



TRPC4AP Rabbit pAb

Catalog No	YP-Ab-18948
Isotype	IgG
Reactivity	
Applications	
Gene Name	
Protein Name	
Immunogen	Synthesized peptide derived from human TRPC4AP
Specificity	This antibody detects endogenous levels of TRPC4AP at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	88kD
Cell Pathway	Cytoplasm, perinuclear region .
Tissue Specificity	
Function	Substrate-recognition component of a DCX (DDB1-CUL4-X-box) E3 ubiquitin-protein ligase complex required for cell cycle control . The DCX(TRPC4AP) complex specifically mediates the polyubiquitination and subsequent degradation of MYC as part of the DesCEND (destruction via C-end degrons) pathway . The DesCEND (destruction via C-end degrons) pathway recognizes a C-degron located at the extreme C terminus of target proteins, leading to their ubiquitination and degradation . The DCX(TRPC4AP) complex specifically recognizes proteins with an arginine at the minus 3 position (R-3 motif) at the C-terminus, such as MYC, leading to their ubiquitination and degradation . Also participates in the activation of NFκB1 in response to ligation of TNFRSF1A, possibly by linking TNFRSF1A to the IKK signalosome (By similarity). Involved in JNK activation via its interaction with TRAF2 (By similarity). Also involved in elevation of endoplasmic reticulum Ca(2+) storage reduction in response



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images