



TRDMT mouse mAb

Catalog No	YP-mAb-09589
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRDMT1 DNMT2
Protein Name	tRNA (cytosine(38)-C(5))-methyltransferase (EC 2.1.1.204) (DNA (cytosine-5)-methyltransferase-like protein 2) (Dnmt2) (DNA methyltransferase homolog HsaIIP) (DNA MTase homolog HsaIIP) (M.HsaIIP) (PuMe
Immunogen	Synthesized peptide derived from human TRDMT AA range: 317-367
Specificity	This antibody detects endogenous levels of human TRDMT
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Observed Band	
Calculated Molecular Weight	43kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Ubiquitous. Higher expression in testis, ovary and thymus and at much lower levels in spleen, prostate, colon, small intestine, and peripheral blood leukocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

