



CASA1 mouse mAb

Catalog No	YP-mAb-09701
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CSN1S1 CASA CSN1
Protein Name	Alpha-S1-casein [Cleaved into: Casoxin-D]
Immunogen	Synthesized peptide derived from human CASA1 AA range: 45-95
Specificity	This antibody detects endogenous levels of human CASA1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	20kD
Cell Pathway	Secreted.
Tissue Specificity	Mammary gland specific. Secreted in milk.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images