



ARRB2 mouse mAb

Catalog No	YP-mAb-17718
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	ARRB2
Protein Name	
Immunogen	Synthetic peptide of human ARRB2
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	arrestin beta 2; ARB2; ARR2; BARR2
Observed Band	Calculated MW: 46 kDa; Observed MW: 46 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Functions in regulating agonist-mediated G-protein coupled receptor (GPCR) signaling by mediating both receptor desensitization and resensitization processes.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images