



c-Jun Mouse mAb

Catalog No	YP-mAb-17706
Isotype	Mouse IgG1
Reactivity	Human, Mouse, Monkey, Rat
Applications	WB, IHC-P, ICC/IF, FC, ELISA
Gene Name	c-Jun
Alternative Names	AP1; AP-1; c-Jun; Jun
Research Field	Epigenetics and Nuclear Signaling
Product Categories	Primary antibody
Host	Mouse
Dilution Ratio	WB: 1/500-1/1000 IHC: 1/100-1/200 IF: 1/50-1/200 FC: 1/50-1/100 ELISA: 1/10000
Molecular Weight	Calculated MW: 43kDa; Observed MW: 43kDa
Clonality	Monoclonal Antibody
Clonality No.	4I7-A9-G2
Immunogen	Purified recombinant fragment of human c-Jun expressed in E. Coli.
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

