



优品生物
UpingBio

网址: www.upingBio.com

邮箱: youpin@youpinbio.com

技术电话: 400-999-8863

Goat Anti Rabbit IgG (Dylight 800)说明书

首先感谢您选择优品生物，我司为您提供优质的免疫学试剂，在您 使用本产品前，请认真阅读说明书，不清楚之处，请联系技术人员解答，这对您顺利完成实 验很有帮助。

Name: Goat Anti Rabbit IgG (Dylight 800)

Catalog No : YPE800103

Applications: IF Flow Cyt WB

Dilution : Suggested starting 1:50-1:1000 dilutions for most fluorescent applications; while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR® Odyssey equipment would be recommended.

Formulation: Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% BSA as stablizer and 50% Glycerol..

Coupling : Dylight 800

Purification: Affinity purified using solid phase Rabbit IgG (H&L) with finally > 95% purity based on SDS-PAGE

Storage Stability: Stable for one year at -20 。 C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

DyLight 350	353/432	Blue
DyLight 405	400/420	Blue
DyLight 488	493/518	Green
Dylight 549	562/576	Yellow
Dylight 594	593/618	Red/Orange
Dylight 649	652/672	Red
Dylight 680	692/712	Near IR
Dylight 800	777/794	Near IR