




AK7 Mouse mAb

Catalog No	YP-mAb-18999
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Protein Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AK 7; FAP75; CFAP75; SPGF27; AK7
Observed Band	87kDa
Calculated Molecular Weight	83kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	This gene encodes a member of the adenylate kinase family of enzymes. The encoded enzyme is a phosphotransferase that catalyzes the reversible phosphorylation of adenine nucleotides. This enzyme plays a role in energy homeostasis of the cell. Alternative splicing results in multiple transcript variants. Mutations in the mouse gene are associated with primary ciliary dyskinesia.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of various lysates using AK7 Mouse mAb (A14600) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25μg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.