



PDIA6 Mouse mAb

Catalog No	YP-mAb-19002
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	PDIA6 ERP5 P5 TXNDC7
Protein Name	Protein disulfide-isomerase A6 (Endoplasmic reticulum protein 5) (ER protein 5) (ERp5) (Protein disulfide isomerase P5) (Thioredoxin domain-containing protein 7)
Immunogen	Synthesized peptide derived from human PDIA6
Specificity	This antibody detects endogenous levels of PDIA6 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	48kD
Cell Pathway	Endoplasmic reticulum lumen . Cell membrane . Melanosome . Identified by mass spectrometry in melanosome fractions from stage I to stage IV (PubMed:12643545). .
Tissue Specificity	Expressed in platelets (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images