



AGAL rabbit pAb

Catalog No	YP-Ab-09486
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	GLA
Protein Name	Alpha-galactosidase A (EC 3.2.1.22) (Alpha-D-galactosidase A) (Alpha-D-galactoside galactohydrolase) (Melibiase) (Agalsidase)
Immunogen	Synthesized peptide derived from human AGAL AA range: 138-188
Specificity	This antibody detects endogenous levels of human AGAL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Calculated Molecular Weight	47kD
Cell Pathway	Lysosome.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

