



XAF1 mouse mAb

Catalog No	YP-mAb-09207
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	XAF1 BIRC4BP XIAPAF1
Protein Name	XIAP-associated factor 1 (BIRC4-binding protein)
Immunogen	Synthesized peptide derived from human XAF1 AA range: 47-97
Specificity	This antibody detects endogenous levels of human XAF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Observed Band	
Calculated Molecular Weight	33kD
Cell Pathway	Cytoplasm. Nucleus. Mitochondrion. Found in the cytoplasm and nucleus of placental syncytiotrophoblasts. Translocates to mitochondria upon TNF-alpha treatment.; [Isoform 1]: Nucleus.; [Isoform 5]: Nucleus.
Tissue Specificity	Widely expressed. Expression is frequently down-regulated in cancer cell lines. Isoform 5 is widely expressed. Expressed in placenta (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

