




PEA3 Monoclonal Antibody

Catalog No	YP-Ab-00082
Isotype	IgG
Reactivity	Human
Applications	WB;ELISA
Gene Name	ETV4
Protein Name	ETS translocation variant 4
Immunogen	Purified recombinant fragment of human PEA3 (aa50-109) expressed in E. Coli.
Specificity	PEA3 Monoclonal Antibody detects endogenous levels of PEA3 protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ETV4; E1AF; PEA3; ETS translocation variant 4; Adenovirus E1A enhancer-binding protein; E1A-F; Polyomavirus enhancer activator 3 homolog; Protein PEA3
Calculated Molecular Weight	54kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in keratinocytes.
Function	function: Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3', similarity: Belongs to the ETS family., similarity: Contains 1 ETS DNA-binding domain.,
Background	function: Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3', similarity: Belongs to the ETS family., similarity: Contains 1 ETS DNA-binding domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis using PEA3 Monoclonal Antibody against K562 cell lysate (1).