



TSPAN2 Mouse mAb

Catalog No	YP-mAb-19022
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	TSPAN2
Protein Name	Tetraspanin-2 (Tspan-2) (Tetraspan NET-3)
Immunogen	Synthesized peptide derived from human TSPAN2
Specificity	This antibody detects endogenous levels of TSPAN2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	24kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	May play a role in signalling in oligodendrocytes in the early stages of their terminal differentiation into myelin-forming glia and may also function in stabilizing the mature sheath.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images