




ACSL4 mouse mAb

Catalog No	YP-mAb-12589
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ACSL4 ACS4 FACL4 LACS4
Protein Name	ACSL4
Immunogen	Synthesized peptide derived from human ACSL4
Specificity	This antibody detects endogenous levels of ACSL4 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.81% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Long-chain-fatty-acid--CoA ligase 4 (EC 6.2.1.3) (Long-chain acyl-CoA synthetase 4) (LACS 4)
Observed Band	75kD
Cell Pathway	Mitochondrion outer membrane ; Single-pass type III membrane protein . Peroxisome membrane ; Single-pass type III membrane protein . Microsome membrane ; Single-pass type III membrane protein . Endoplasmic reticulum membrane ; Single-pass type III membrane protein . Cell membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Regulation of Fatty Acid Metabolism and Inhibition of Colorectal Cancer Progression by Erchen Decoction Evidence-based Complementary and Alternative

Medicine Linghong Liao, Fei Zhang, Zewei Zhuo, Chengbao Huang, Xiaofang Zhang, Ruifang Liu, Bizhen Gao, Shanshan Ding WB Mouse colorectal tissue

