



DUS8 rabbit pAb

Catalog No	YP-Ab-09999
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	DUSP8 C11orf81 VH5
Protein Name	Dual specificity protein phosphatase 8 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase hVH-5)
Immunogen	Synthesized peptide derived from human DUS8 AA range: 98-148
Specificity	This antibody detects endogenous levels of human DUS8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Calculated Molecular Weight	69kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Abundant in brain, heart and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

