



FBXL15 Rabbit pAb

Catalog No	YP-Ab-19070
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	FBXL15 FBXO37
Protein Name	F-box/LRR-repeat protein 15 (F-box only protein 37)
Immunogen	Synthesized peptide derived from human FBXL15
Specificity	This antibody detects endogenous levels of FBXL15 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	33kD
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	Substrate recognition component of a SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of SMURF1, thereby acting as a positive regulator of the BMP signaling pathway. Required for dorsal/ventral pattern formation and bone mass maintenance. Also mediates ubiquitination of SMURF2 and WWP2.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images