



GABP-α Mouse mAb

Catalog No	YP-mAb-00982
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GABPA
Protein Name	GA-binding protein alpha chain
Immunogen	Purified recombinant fragment of human GABP-α (aa120-190) expressed in E. Coli.
Specificity	GABP-α Monoclonal Antibody detects endogenous levels of GABP-α protein.
Formulation	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GABPA; E4TF1A; GA-binding protein alpha chain; GABP subunit alpha; Nuclear respiratory factor 2 subunit alpha; Transcription factor E4TF1-60
Observed Band	51kD
Cell Pathway	Nucleus.
Tissue Specificity	Brain,Hepatoma,
Function	function:Transcription factor capable of interacting with purine rich repeats (GA repeats). Necessary for the expression of the Adenovirus E4 gene.,similarity:Belongs to the ETS family.,similarity:Contains 1 ETS DNA-binding domain.,similarity:Contains 1 PNT (pointed) domain.,subunit:Heterotetramer of two alpha and two beta subunits.,
Background	This gene encodes one of three GA-binding protein transcription factor subunits which functions as a DNA-binding subunit. Since this subunit shares identity with a subunit encoding the nuclear respiratory factor 2 gene, it is likely involved in activation of cytochrome oxidase expression and nuclear control of mitochondrial function. This subunit also shares identity with a subunit constituting the transcription factor E4TF1, responsible for expression of the adenovirus E4 gene. Because of its chromosomal localization and ability to form heterodimers with other polypeptides, this gene may play a role in the Down Syndrome phenotype. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Oct 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis using GABP- α Monoclonal Antibody against HeLa (1), A549 (2), MCF-7 (3), NIH/3T3 (4) and SMMC-7721 (5) cell lysate.