



EVX1 mouse mAb

Catalog No	YP-mAb-09277
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EVX1
Protein Name	Homeobox even-skipped homolog protein 1 (EVX-1)
Immunogen	Synthesized peptide derived from human EVX1 AA range: 312-362
Specificity	This antibody detects endogenous levels of human EVX1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	45kD
Cell Pathway	Nucleus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

