



mon1b Rabbit pAb

Catalog No	YP-Ab-19100
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	MON1B HSRG1 KIAA0872 SAND2
Protein Name	Vacuolar fusion protein MON1 homolog B (HSV-1 stimulation-related gene 1 protein) (HSV-I stimulating-related protein)
Immunogen	Synthesized peptide derived from human MON1B
Specificity	This antibody detects endogenous levels of MON1B at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	60kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images