



FBSP1 mouse mAb

Catalog No	YP-mAb-09210
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FBXO45 FBX45
Protein Name	F-box/SPRY domain-containing protein 1 (F-box only protein 45) (hFbxo45)
Immunogen	Synthesized peptide derived from human FBSP1 AA range: 120-170
Specificity	This antibody detects endogenous levels of human FBSP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	31kD
Observed Band	
Cell Pathway	Cell junction, synapse, postsynaptic cell membrane . Cell junction, synapse, presynaptic cell membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

