



BIEA rabbit pAb

Catalog No	YP-Ab-10020
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BLVRA BLVR BVR
Protein Name	Biliverdin reductase A (BVR A) (EC 1.3.1.24) (Biliverdin-IX alpha-reductase)
Immunogen	Synthesized peptide derived from human BIEA AA range: 84-134
Specificity	This antibody detects endogenous levels of human BIEA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	33kD
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

