



FA83A mouse mAb

Catalog No	YP-mAb-09528
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FAM83A TSGP
Protein Name	Protein FAM83A (Tumor antigen BJ-TSA-9) (Tumor-specific gene expressed in prostate protein)
Immunogen	Synthesized peptide derived from human FA83A AA range: 374-424
Specificity	This antibody detects endogenous levels of human FA83A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	48kD
Observed Band	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

