



YN4453

TM155 mouse mAb

Catalog No	YP-mAb-09293
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TMEM155
Protein Name	Protein TMEM155
Immunogen	Synthesized peptide derived from human TM155 AA range: 31-81
Specificity	This antibody detects endogenous levels of human TM155
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	14kD
Observed Band	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

