



ZMPSTE24 Rabbit pAb

Catalog No	YP-Ab-19142
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	ZMPSTE24 FACE1 STE24
Protein Name	CAAX prenyl protease 1 homolog (Farnesylated proteins-converting enzyme 1) (FACE-1) (Prenyl protein-specific endoprotease 1) (Zinc metalloproteinase Ste24 homolog)
Immunogen	Synthesized peptide derived from human ZMPSTE24
Specificity	This antibody detects endogenous levels of ZMPSTE24 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	52kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Nucleus inner membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed. High levels in kidney, prostate, testis and ovary.
Function	Proteolytically removes the C-terminal three residues of farnesylated proteins. Acts on lamin A/C.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

