



NA14 Rabbit pAb

Catalog No	YP-Ab-19157
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	SSNA1 NA14
Protein Name	Sjogren syndrome nuclear autoantigen 1 (Nuclear autoantigen of 14 kDa)
Immunogen	Synthesized peptide derived from human NA14
Specificity	This antibody detects endogenous levels of NA14 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	13kD
Cell Pathway	Nucleus. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

