



GFI1B mouse mAb

Catalog No	YP-mAb-12423
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GFI1B
Protein Name	GFI1B
Immunogen	Synthesized peptide derived from human GFI1B
Specificity	This antibody detects endogenous levels of Human GFI1B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	37kD
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

