



Vstm2L Mouse mAb

Catalog No	YP-mAb-19174
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	VSTM2L C20orf102
Protein Name	V-set and transmembrane domain-containing protein 2-like protein
Immunogen	Synthesized peptide derived from human VSTM2L
Specificity	This antibody detects endogenous levels of VSTM2L at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	22kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images