



## FBXW4 Mouse mAb

Catalog No	YP-mAb-19199
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FBXW4 FBW4 SHFM3
Protein Name	F-box/WD repeat-containing protein 4 (Dactylin) (F-box and WD-40 domain-containing protein 4)
Immunogen	Synthesized peptide derived from human FBXW4
Specificity	This antibody detects endogenous levels of FBXW4 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	45kD
Cell Pathway	Expressed in brain, kidney, lung and liver.
Tissue Specificity	
Function	Probably recognizes and binds to some phosphorylated proteins and promotes their ubiquitination and degradation. Likely to be involved in key signaling pathways crucial for normal limb development. May participate in Wnt signaling.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images