



WDR92 Rabbit pAb

Catalog No	YP-Ab-19210
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WDR92
Protein Name	WD repeat-containing protein 92 (WD repeat-containing protein Monad)
Immunogen	Synthesized peptide derived from human WDR92
Specificity	This antibody detects endogenous levels of WDR92 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	39kD
Cell Pathway	Dynein axonemal particle .
Tissue Specificity	Widely expressed with the highest expression in testis.
Function	Key assembly factor specifically required for the stability of axonemal dynein heavy chains in cytoplasm.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images