




MET14 mouse mAb

Catalog No	YP-mAb-12013
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	METTL14 KIAA1627
Protein Name	MET14
Immunogen	Synthesized peptide derived from human MET14 AA range: 360-410
Specificity	This antibody detects endogenous levels of MET14 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	50kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:Probable methyltransferase.,similarity:Belongs to the MT-A70-like family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night