



KRA31 rabbit pAb

Catalog No	YP-Ab-09729
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KRTAP3-1 KAP3.1 KRTAP3.1
Protein Name	Keratin-associated protein 3-1 (High sulfur keratin-associated protein 3.1) (Keratin-associated protein 3.1)
Immunogen	Synthesized peptide derived from human KRA31 AA range: 11-61
Specificity	This antibody detects endogenous levels of human KRA31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	11kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

